

# Mitochondrial respiration and ATP production <sup>[1, 2]</sup> on homogenate of red muscle of Sea bass

## Products:

- **MiR05 buffer** (without EGTA and MgCl<sub>2</sub>) - aliquots (20\*50 mL) at -20°C

	Sigma ref	Store at	Final conc.	FW	Add to 1L
HPLC H <sub>2</sub> O					500 mL
K-lactobionate	153516-100G	RT	60mM	358.3	120mL of 0.5M*
Taurine	T0625-25G	RT	20mM	125.1	2.502g
KH <sub>2</sub> PO <sub>4</sub>	P9791-1KG	RT	10mM	136.1	1.361g
HEPES	H7523-250G	RT	20mM	238.3	4.77g
Sucrose	S9378-1KG	RT	110mM	342.3	37.65g
BSA <sup>1</sup>	A6003-25G	4 °C	1g/L		1g
<sup>2</sup> KOH (Ca <sup>2+</sup> free)	P5958-250G	RT	5N	56.11	pH 7.2
HPLC H <sub>2</sub> O					Qsp 1L

\*K-lactobionate 0.5M - Freshly prepared

	Vol final = 240 mL	Vol final = 120 mL
Lactobionic acid	43 g	21.50 g
HPLC H <sub>2</sub> O	160 mL	80 mL
KOH (Ca <sup>2+</sup> free)	pH 7.0 at RT	pH 7.0 at RT
HPLC H <sub>2</sub> O	Qsp 240 mL	Qsp 120 mL

- **Adenylate kinase inhibitor (Ap5A)** - aliquots of 100µL at -20°C (07/04/2019 - Lot SLBT4515)

	Sigma ref	Store at	Final conc.	FW	1950µL dH <sub>2</sub> O
Ap5A	D4022-50mg	-20°C	25mM	1026.28	50mg

- **Pyruvate** - Freshly prepared

	Sigma ref	Store at	Final conc.	FW	200µL H <sub>2</sub> O δ
Sodium pyruvate	P2256-25 g	4°C	2M	110.0	44 mg

- **Malate** - aliquots (10\*500 µL) at -20°C

	Sigma ref	Store at	Final conc.	FW	5mL H <sub>2</sub> O δ
Malic acid	M1000-100 g	RT	400mM	134.1	268.2mg
HPLC H <sub>2</sub> O					3mL
KOH (Ca <sup>2+</sup> free)		RT	5N	56.11	pH 7.1 at RT (≈ 700 µL)
HPLC H <sub>2</sub> O					Qsp 5mL

<sup>1</sup> Preliminary dissolved in heated up H<sub>2</sub>O δ

<sup>2</sup> No longer sold by Sigma, found under same ref at [Honeywell Speciality Chemicals](#), Germany

- Succinate disodium salt, hexahydrate - aliquots (10\*500  $\mu$ L) at -20°C

	Sigma ref	Store at	Final conc.	FW	5mL
Succinate	S2378-100 g	RT	1M	270.1	1.3505g
HPLC H <sub>2</sub> O					3mL
TrisHCl	T3253-250 g	RT	1N	157.60	pH 7
HPLC H <sub>2</sub> O					Qsp 5mL

- Magnesium Green 5N hexapotassium salt- aliquots (20\*50 $\mu$ L; in opaque tubes) at -20°C

	Life technologies	Store at	Final conc.	FW	1mL HPLC H <sub>2</sub> O
MgG 5K <sup>+</sup>	M3733-1mg	-20°C	1.1mM	915.9004	1mg

- Magnesium chloride (MgCl<sub>2</sub>) - Room temperature

	Sigma ref	Store at	Final conc.	FW	10mL HPLC H <sub>2</sub> O
MgCl <sub>2</sub> •6H <sub>2</sub> O	M2670-500g	RT	1M	203.30	2.0g

- Dilute to 200mM by adding 2mL of 1M stock to 8mL of HPLC H<sub>2</sub>O.

- EDTA - Room temperature

	Sigma ref	Store at	Final conc.	FW	10mL HPLC H <sub>2</sub> O
EDTA	E9884-100g	RT	200mM	292.24	585mg
KOH	P5958-250g	RT	5N	56.11	pH 6
Heat up solution					
KOH	P5958-250g	RT	5N	56.11	pH 7.1

- Dilute to 1mM by adding 50 $\mu$ L of 200mM stock to 9950 $\mu$ L of HPLC H<sub>2</sub>O.

- EGTA - Room temperature

	Sigma ref	Store at	Final conc.	FW	10mL HPLC H <sub>2</sub> O
EGTA	E4378-100mg	RT	200mM	380.35	761mg
KOH	P5958-250g	RT	5N	56.11	pH 7.1
Heat up solution					
KOH	P5958-250g	RT	5N	56.11	pH 7.1

- ADP K<sup>+</sup> 250 mM - aliquots (40\*200 $\mu$ L) at -80°C (NO freeze-thaw cycle)

	Merck ref	Store at	Final conc.	FW	8mL HPLC H <sub>2</sub> O
ADP	117105-1G	-20°C	250 mM	501.3	1g
HPLC H <sub>2</sub> O					4mL
KOH	P5958-250G	RT	5 N	56.11	pH 6.3-6.9 at RT (~750 $\mu$ L)
HPLC H <sub>2</sub> O					Qsp 8 mL

- **Carboxylatractyloside potassium salt from *Xanthium sibiricum*** - aliquots (11\*200 $\mu$ L) at -20°C

	Sigma ref	Store at	Final conc.	FW	2.1mL HPLC H <sub>2</sub> O
<b>cATR</b>	216201-2mg	-20°C	1mM	939.1	2mg

- **Antimycin A** (from *Streptomyces* sp.) - aliquots of 200 $\mu$ L at -20°C.

	Sigma ref	Store at	Final conc.	FW	2mL Eth 99%
<b>Antimycin A</b>	A8674-25 mg	-20 °C	5mM	540.0	5.4mg

- **Ascorbate sodium salt** (Sodium L-ascorbate) - aliquots (20\*500 $\mu$ L) in opaque tubes in -20°C - use one a day

	Sigma ref	Store at	Final conc.	FW	10mL HPLC H <sub>2</sub> O
<b>Ascorbate</b>	A4034-100 g	RT	800mM	198.11	1.584g

- **N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride**- aliquots (10\*120 $\mu$ L) in opaque tubes at -20°C - use one a day

	Sigma ref	Store at	Final conc.	FW	1.2mL HPLC H <sub>2</sub> O
<b>TMPD</b>	T3134-5 g	RT	200mM	237.2	56.9mg

#### Cleaning:

- 3 times distilled water
- 2 times Eth 70%
- 5 min Eth 70%
- 10 min Eth 100%
- 5 min Eth 100% ----- or: Overnight Eth 70%
- 2 times Eth 70%
- 3 times distilled water

#### References:

1. Salin, K., et al., *Variation in metabolic rate among individuals is related to tissue-specific differences in mitochondrial leak respiration*. *Physiological and Biochemical Zoology*, 2016. **89**(6): p. 511-523.
2. Salin, K., et al., *Simultaneous measurement of mitochondrial respiration and ATP production in tissue homogenates and calculation of effective P/O ratios*. *Physiological Reports*, 2016. **4**(20): p. e13007.

**Methods:**

-2.1mL of **red muscle** (1 mg/mL) homogenate in MiRO5 buffer (without EGTA and  $\text{MgCl}_2$ ).

-----  
 -Add 8 $\mu\text{L}$  of **Ap5A** 25mM (Cf = 100 $\mu\text{M}$ ) with a pipette

-----Put stoppers-----

Chamber \_ Oroboros \_  $\rightarrow$  Add oxygen the medium (550  $\mu\text{M}$  <  $[\text{O}_2]_{13^\circ\text{C}}$  < 650  $\mu\text{M}$ )

-----OR-----

Chamber \_ Oroboros \_  $\rightarrow$  Add oxygen the medium (550  $\mu\text{M}$  <  $[\text{O}_2]_{16^\circ\text{C}}$  < 650  $\mu\text{M}$ )

-----OR-----

Chamber \_ Oroboros \_  $\rightarrow$  NO oxygen the medium ( $[\text{O}_2]_{16^\circ\text{C}}$  = 276  $\mu\text{M}$ )

-----OR-----

Chamber \_ Oroboros \_  $\rightarrow$  NO oxygen the medium ( $[\text{O}_2]_{13^\circ\text{C}}$  = 295  $\mu\text{M}$ )

-----Close chamber-----

- Add 5 $\mu\text{L}$  of **Pyruvate** 2M (Cf = 5mM)

- Add 2.5 $\mu\text{L}$  of **Malate** 400mM (Cf = 0.5mM)

-----  
 - Add 4 $\mu\text{L}$  of **MgGreen** 1.1mM (Cf = 2.2 $\mu\text{M}$ )

-----Reading for 30sec-----

-Add 1 $\mu\text{L}$  of **EGTA** 200mM (Cf = 0.1mM)

-----Reading for 30sec -----

-Add 10 $\mu\text{L}$  of **EDTA** 1mM (Cf = 5 $\mu\text{M}$ )

-----Reading for 30sec -----

REPEAT PREVIOUS STEP 2 MORE TIMES [EDTA (Cf = 15 $\mu\text{M}$ )]

-----  
 -Add 1 $\mu\text{L}$  of **MgCl<sub>2</sub>** 200mM (Cf = 0.1mM)

-----Reading for 30sec -----

REPEAT PREVIOUS STEP 9 MORE TIMES [MgCl<sub>2</sub> (Cf = 1mM)]

-----Reading for 30sec -----

-Add 20 $\mu\text{L}$  of **Succinate** 1M (Cf = 10 $\mu\text{M}$ )

-----Reading for 5 minutes-----

-Add 16 $\mu\text{L}$  of **ADP** 250mM (Cf = 2mM)

----- Reading until stabilization -----

-Add 8 $\mu\text{L}$  of **cATR** 1mM (Cf = 4 $\mu\text{M}$ )

-----Reading until stabilization-----

- Add 1 $\mu\text{L}$  of **Antimycin A** 5mM (Cf = 2.5 $\mu\text{M}$ )

-----Reading until stabilization-----

- Add 20 $\mu\text{L}$  of **Ascorbate** 800mM (Cf = 8mM)

-----Reading until stabilization-----

- Add 5 $\mu\text{L}$  of **TMPD** 200mM (Cf = 500 $\mu\text{M}$ )

-----Reading 5 minutes-----